CageControl

Control waveplates inside tomography cages

Generated by Doxygen 1.8.13

Contents

1	Bug	List													1
2	Nam	espace	Index												3
	2.1	Names	space List				 		 3						
3	Hier	archica	l Index												5
	3.1	Class	Hierarchy				 		 5						
4	Clas	s Index	:												7
	4.1	Class	List				 		 7						
5	File	Index													9
	5.1	File Lis	st				 		 9						
6	Nam	nespace	Docume	ntation											11
	6.1	helper	Namespa	ce Refer	rence		 		 11						
		6.1.1	Detailed	Descrip	tion		 		 11						
		6.1.2	Function	Docum	entatio	on .	 		 11						
			6.1.2.1	error()			 		 11						
			6.1.2.2	info()			 		 12						
			6.1.2.3	messa	ige()		 		 12						
			6.1.2.4	warnin	ıg() .		 		 12						

ii CONTENTS

7	Clas	s Docu	mentation		13
	7.1	cageco	ontrol Class	s Reference	13
		7.1.1	Member F	Function Documentation	16
			7.1.1.1	initconnections()	17
			7.1.1.2	LoadConfig()	17
			7.1.1.3	motorGB()	17
			7.1.1.4	movemotor()	17
			7.1.1.5	SaveConfig()	18
			7.1.1.6	slot_changeoffsetusage	18
			7.1.1.7	slot_changeWPangles	18
			7.1.1.8	slot_moveHV	18
			7.1.1.9	slot_moveLR	19
			7.1.1.10	slot_movemotors	19
			7.1.1.11	slot_movePM	19
			7.1.1.12	updatesettings	20
			7.1.1.13	updatesettingsint	20
			7.1.1.14	updatestatus()	20
		7.1.2	Member [Data Documentation	20
			7.1.2.1	useoffset	21
	7.2	Motor	Class Refe	rence	21
		7.2.1	Detailed [Description	23
		7.2.2	Member F	Function Documentation	23
			7.2.2.1	command_microstep	23
			7.2.2.2	command_moveboth	24
			7.2.2.3	command_singlestep	24
			7.2.2.4	command_step	24
			7.2.2.5	handleError	25
			7.2.2.6	isopen	25
			7.2.2.7	motorstatusmessage	25
			7.2.2.8	sensordata()	25

CONTENTS

		8.6.1	Detailed	Description	41
	8.6	•		erence	40
	8.5			Reference	39
	8.4			rence	38
	8.3	helper.		erence	37
			8.2.2.9	UNUSED	37
			8.2.2.8	RADTODEG	37
			8.2.2.7	PI	37
			8.2.2.6	EPS	36
			8.2.2.5	DEGTORAD	36
			8.2.2.4	DEBUGWARNING	36
			8.2.2.3	DEBUGINFO	36
			8.2.2.1	DEBUGERROR	36
		0.2.2	8.2.2.1	DEBUG	36 36
		8.2.1		Description	35
	8.2	8.2.1		erence	34
	0.0	8.1.1		Description	34
	8.1			Prence	33
8			entation		33
_					
			7.3.3.8	processCommands	31
			7.3.3.7	MovePM	31
			7.3.3.6	MoveLR	30
			7.3.3.5	MoveHV	30
			7.3.3.4	Move	30
			7.3.3.3	changeWPangles	30
			7.3.3.1	changeoffsetusage	29
		1.3.3	7.3.3.1	bind	29 29
		7.3.3	7.3.2.1 Member	Function Documentation	29 29
		7.3.2		tor & Destructor Documentation	29 29
		7.3.1		Description	28
	7.3			s Reference	27
	7.0	LIDD"	7.2.3.1	hometimer	27
		7.2.3		Data Documentation	26
			7.2.2.11	write	26
			7.2.2.10	stop	26
			7.2.2.9	showStatusMessage	26

Bug List

File debug.h

Printing to console does not work on Windows. Workaround: Redirect stderr to stdout and redirect stdout to a file.

File defines.h

There are no known bugs.

Namespace helper

There are no known bugs.

Class Motor

There are no known bugs.

Class UDPlistener

There are no known bugs

File version.h

There are no known bugs.

2 Bug List

Namespace Index

2.1	Namespace	Liet
4.1	Ivallicapace	LIS

Here is a list of all documented namespaces with brief descrip-	olions:

helper																
	Small functions to display messages			 												1

4 Namespace Index

Hierarchical Index

3.1 Class Hierarchy

This inheritance list is sorted roughly, but not completely, alphabetically:

QMainWindow																		
cagecontrol	 		 								 							 13
QObject																		
Motor	 		 								 							 21
LIDPlistaner																		27

6 Hierarchical Index

Class Index

4.1 Class List

Here are the classes, structs, unions and interfaces with brief descriptions:

cagecontrol	 	 13
Motor		
Operates the PCB-motor	 	 21
UDPlistener		
Used to control dinspect with UDP packages	 	 27

8 Class Index

File Index

5.1 File List

Here is a list of all documented files with brief descriptions:

cagecontrol.h	??
debug.h	
Debug macros	33
defines.h	
Various compile-time definitions	34
helper.h	37
motor.h	
udplistener.h	39
version.h	
This file contains information about the code version	40

10 File Index

Namespace Documentation

6.1 helper Namespace Reference

contains small functions to display messages

Functions

• void message (QString msg)

message displays a message box

• void error (QString msg)

error displays an error-messagebox and writes a debug_error message to stdout

• void warning (QString msg)

warning displays warning-messagebox and writes a debug_warning message to stdout

• void info (QString msg)

info displays an info-messagebox and writes a debug_info message to stdout

6.1.1 Detailed Description

contains small functions to display messages

Bug There are no known bugs.

6.1.2 Function Documentation

6.1.2.1 error()

error displays an error-messagebox and writes a debug_error message to stdout

Parameters

msg the message to be displayed

6.1.2.2 info()

```
void helper::info (
          QString msg )
```

info displays an info-messagebox and writes a debug_info message to stdout

Parameters

msg the message to be displayed

6.1.2.3 message()

message displays a message box

Parameters

msg the message to be displayed

6.1.2.4 warning()

warning displays warning-messagebox and writes a debug_warning message to stdout


Parameters

msg the message to be displayed


Class Documentation

7.1 cagecontrol Class Reference

Inheritance diagram for cagecontrol:



Collaboration diagram for cagecontrol:



Public Slots

void slot_changeWPangles (QVector< double > angles)

slot_changeWPangles sets offsetangles for all waveplates

· void slot changeoffsetusage (bool uo in)

slot_changeoffsetusage changes the usage of the waveplate offset

void slot_moveHV (QString color)

moveHV

void slot movePM (QString color)

movePM

void slot_moveLR (QString color)

moveLF

• void slot_movemotors (QString color, double HWPang, double QWPang)

movemotors

Public Member Functions

cagecontrol (QWidget *parent=nullptr)

Private Slots

· void updatesettings (double d)

updatesettings fills variables with data from GUI

void updatesettingsint (int i)

updatesettingsint wrapper, just calls

void updateUI ()

updateUI updates UI with supposedly new numbers (loaded from conf file, e.g.)

Private Member Functions

void setupUI (QGridLayout *layout)

Puts together the GUI.

• void openmotors ()

Opens serial connections to the PCB motor controlllers.

void updatestatus (QString msg)

updatestatus writes message to statusbar and to a logfile

void LoadConfig ()

LoadConfig loads config from a file.

• void SaveConfig ()

SaveConfig stores config to a file.

void motorGB (QGroupBox *gb, QString id)

motorGB fills an empty QGroupBox with motor controls

• void initconnections ()

initconnections connects Qt Signals to slots

void movemotor (QString motor, double HWPang, double QWPang)

movemotor moves both motors of a cage to certain angles

void moveredHV ()

moveredHV moves red cage to HV basis

void moveredPM ()

```
moveredPM moves red cage to PM basis
· void moveredLR ()
     moveredLR moves red cage to RL basis
· void moveredANG ()
     moveredANG moves red cage to the angles set in the GUI

    void movebrownHV ()

     movebrownHV moves brown cage to HV basis

    void movebrownPM ()

     movebrownPM moves brown cage to PM basis

    void movebrownLR ()

     movebrownLR moves brown cage to RL basis

    void movebrownANG ()

     movebrownANG moves brown cage to the angles set in the GUI

    void movegreenHV ()

     movegreenHV moves green cage to HV basis

    void movegreenPM ()

     movegreenPM moves green cage to PM basis

    void movegreenLR ()

     movegreenLR moves green cage to RL basis
• void movegreenANG ()
     movegreenANG moves green cage to the angles set in the GUI

    void moveblueHV ()

     moveblueHV moves blue cage to HV basis

    void movebluePM ()

     movebluePM moves blue cage to PM basis

    void moveblueLR ()

     moveblueLR moves blue cage to RL basis
• void moveblueANG ()
     moveblueANG moves white cage to the angles set in the GUI

    void movewhiteHV ()

     movewhiteHV moves white cage to HV basis

    void movewhitePM ()

     movewhitePM moves white cage to PM basis

    void movewhiteLR ()

     movewhiteLR moves white cage to RL basis

    void movewhiteANG ()

     movewhiteANG moves black cage to the angles set in the GUI

    void moveblackHV ()

     moveblackHV moves black cage to HV basis

    void moveblackPM ()

     moveblackPM moves black cage to PM basis
· void moveblackLR ()
     moveblackLR moves black cage to RL basis

    void moveblackANG ()

     moveblackANG
· void moveallhv ()
     moveallhy
· void moveallpm ()
     moveallpm
· void movealllr ()
     movealllr

    void moveallarb ()
```

moveallarb

Private Attributes

int udpport

Hold the UDP port to listen to for commandds.

· bool pauseupdating

Keep updateUI and updatesettings from interfering with each other.

· bool useoffset

If true, the angles in the settings-tab will be used as '0'.

QSettings * settings

A QSettings object, used to store settings in a config file.

• UDPlistener * udplistener

Listens to a UDP port, aquiires & checks commands send to it.

QTabWidget * tabs

GUI tab widget.

QWidget * settingstab

GUI tab containing settings.

QWidget * motorstab

GUI tab containing motor controls.

• QStatusBar * status

Status bar.

QVector< QString > comports

Vector containing available serial ports names ports.

QVector< bool > invert

True: invert predefined bases (H/V -> V/H, P/M->M/P, L/R->R/L)

QVector< Motor * > motors

List of serial connections to the cages.

• QVector< QString > motorName

List of colorcodes of the cages.

QVector< QDoubleSpinBox > HWP0sp

List of QSpinBoxes to set the '0' of the HWPs.

 $\bullet \quad \mathsf{QVector}{<\mathsf{QDoubleSpinBox}} > \mathsf{QWP0sp}$

List of QSpinBoxes to set the '0' of the QWPs.

QVector< int > HWPmnum

Motornumber of controller the HWP is connected to.

QVector< int > QWPmnum

Motornumber of controller the QWP is connected to.

• QVector< double > HWP0

'0' of HWPs

• QVector< double > QWP0

'0' of QWPs

• QVector< double > HWPcust

custum set angle to rotate HWP to

• QVector< double > QWPcust

custom set angle to rotate QWP to

QVector< QGroupBox * > uiMotorGroupBoxes

List of Groupboxes containing cage controls.

7.1.1 Member Function Documentation

7.1.1.1 initconnections()

```
void cagecontrol::initconnections ( ) [private]
```

initconnections connects Qt Signals to slots

Defines what happens when a button is clicked, a number is changet, et cetera

7.1.1.2 LoadConfig()

```
void cagecontrol::LoadConfig ( ) [private]
```

LoadConfig loads config from a file.

The dialog is set up with values already stored in the QSettings object. If a specific quantity does not exist there, it is set to a standard value.

7.1.1.3 motorGB()

```
void cagecontrol::motorGB (  \mbox{QGroupBox} * gb, \\ \mbox{QString } id \mbox{)} \mbox{ [private]}
```

motorGB fills an empty QGroupBox with motor controls

Parameters

gb	empty QGroupBox
id	colorcode of the cage

7.1.1.4 movemotor()

movemotor moves both motors of a cage to certain angles

Parameters

motor	colorcode of the cage
HWPang	angle of the HWP in degrees
QWPang	angle of the QWP in degrees

7.1.1.5 SaveConfig()

```
void cagecontrol::SaveConfig ( ) [private]
```

SaveConfig stores config to a file.

The QSettings object is updated with the values received from the dialog and saved immediately.

7.1.1.6 slot_changeoffsetusage

```
void cagecontrol::slot_changeoffsetusage ( bool\ uo\_in\ ) \quad [slot]
```

slot_changeoffsetusage changes the usage of the waveplate offset

Parameters

```
useoffset | true if waveplate offset is to be used
```

The term 'offset' refers to the waveplate angles specified in the

See also

settingstab. E.g.: 'H' of HWP specified in he settingstab is 50° and one wants to rotate the waveplate to H+10°. If (useoffset==true), one needs to rotate the motor to 10° . If (useoffset==false), one needs to rotate the motor to 60° .

7.1.1.7 slot_changeWPangles

```
void cagecontrol::slot_changeWPangles ( {\tt QVector} < {\tt double} \ > {\tt angles} \ ) \quad [{\tt slot}]
```

slot_changeWPangles sets offsetangles for all waveplates

Parameters

```
angles vector containing all angles. Ordering: HWP0,HWP1,...,HWPn,QWP0,QWP1,...,QWPn
```

7.1.1.8 slot_moveHV

moveHV

Parameters

color colorcode of the cage, or 'all'

7.1.1.9 slot_moveLR

moveLR

Parameters

color colorcode of the cage, or 'all'

7.1.1.10 slot_movemotors

movemotors

Parameters

color colorcode of the cage, or 'all'

7.1.1.11 slot_movePM

movePM

Parameters

color colorcode of the cage, or 'all'

7.1.1.12 updatesettings

```
void cagecontrol::updatesettings ( \label{eq:double} \mbox{double } \mbox{$d$ } \mbox{$)$} \ \ [\mbox{private}], \ [\mbox{slot}]
```

updatesettings fills variables with data from GUI

Parameters

d unused

7.1.1.13 updatesettingsint

```
void cagecontrol::updatesettingsint ( \quad \text{ int } i \text{ ) [private], [slot]}
```

updatesettingsint wrapper, just calls

See also

updatesettings(double d)

Parameters

i unused

7.1.1.14 updatestatus()

updatestatus writes message to statusbar and to a logfile

Parameters

msg Message to write

7.1.2 Member Data Documentation

7.2 Motor Class Reference 21

7.1.2.1 useoffset

```
bool cagecontrol::useoffset [private]
```

If true, the angles in the settings-tab will be used as '0'.

See also

slot_changeoffsetusage

The documentation for this class was generated from the following files:


- · cagecontrol.h
- · cagecontrol.cpp
- cagecontrol_motors.cpp
- cagecontrol_ui.cpp

7.2 Motor Class Reference


The Motor class operates the PCB-motor.

```
#include <motor.h>
```

Inheritance diagram for Motor:



Collaboration diagram for Motor:



Public Slots

void open (QString port)

open establishes a connection over a serial port

· void close ()

close closes the serialport connection

· void read ()

read reads from the serial port

void write (const QByteArray &data)

write writes to the serialport

void handleError (QSerialPort::SerialPortError error)

handleError prints an error message of the serialport connection and closes the connection

void showStatusMessage (const QString &message)

showStatusMessage fills the label in the GUI with text

• bool isopen ()

isopen returns the state of the serial connection

void command_park ()

command_park moves the motor to the mechanical stop

void command home ()

command_home sends commands to position at the mechanical stop and afterwards go to the offset starting position, but in an inaccurate way

void command info ()

command info sends the command to request the PCBMotor information

void command help ()

command_help sends the command to print the PCBMotor help

• void command_frequency_sweep ()

command_frequency_sweep sends the PCBMotor command for a frequency sweep

void command_singlestep (QString dirstring)

command_singlestep moves the motor a single step in a direction specified by dirstring

void command_step (uint16_t numsteps, QString dirstring)

command_step moves the motor numstep steps in a direction specified by dirstring

void command microstep (uint16 t nummsteps, QString dirstring)

command microstep aplies nummsteps micropulses to the motor

void stop (bool stop)

stop Tries to stop movenents if possible

void command_moveboth (double ang1, double ang2)

command moveboth moves both motors connected to the controller

Signals

void motorstatusmessage (const QString &message)

motorstatusmessage emitted when the status of the serial connection changes, with a string indicating the actual state.

void ConnectionClosed ()

emitted when serial connection is closed

Public Member Functions

• Motor ()

Motor the contructor initializes variables and establishes the serial connection.

• bool sensordata ()

sensordata returns the current PCBMotor optical encoder wheel sendor state

7.2 Motor Class Reference 23

Public Attributes

QString publicmotorstatusmessage

A string containing the current state of the serial connection.

QSerialPort * serial

Qt serial connection interface.

Private Member Functions

void moveboth ()

command_moveboth moves both motors connected to the controller

Private Attributes

QTimer hometimer

Used to iterate through the steps of 'go to the starting position' - but in an inaccurate way.

· QTimer bothtimer

Used to iterate through the steps of moving two motors of one controller.

· int movebothstep

Controls logic flow when two motors are to be moved consecutively.

· bool serialconnectionok

False if opening the serial connection failed.

• uint16 t motor1steps

number of stept the 1st motor is to be moved

• uint16_t motor2steps

number of stept the 2nd motor is to be moved

7.2.1 Detailed Description

The Motor class operates the PCB-motor.

Bug There are no known bugs.

The PCBMotor is controllable by sending ASCII commands over a serial connection. This class establishes such a connection and controls the movements of the motor.

7.2.2 Member Function Documentation

7.2.2.1 command_microstep

command_microstep aplies nummsteps micropulses to the motor

Parameters

nummsteps	number of micropulses to apply
dirstring	string containing the desired direction

dirstring may either be "bw" of "fw" for backward/forward movement.

7.2.2.2 command_moveboth

command_moveboth moves both motors connected to the controller

Parameters

ang1	angle motor 1 is to be moved to
ang2angle	motor 2 is to be moved to

7.2.2.3 command_singlestep

command_singlestep moves the motor a single step in a direction specified by dirstring

Parameters

dirstrina	a string containing the desired movement direction
acag	

Dirstring may either be "bw" of "fw" for backward/forward movement.

7.2.2.4 command_step

command_step moves the motor numstep steps in a direction specified by dirstring

Parameters

numsteps	number of steps to go
dirstring	direction to go

7.2 Motor Class Reference 25

Dirstring may either be "bw" of "fw" for backward/forward movement.

7.2.2.5 handleError

handleError prints an error message of the serialport connection and closes the connection

Parameters

error

7.2.2.6 isopen

```
bool Motor::isopen ( ) [slot]
```

isopen returns the state of the serial connection

Returns

true if serial connection was established successfully, false otherwise

7.2.2.7 motorstatusmessage

motorstatusmessage emitted when the status of the serial connection changes, with a string indicating the actual state.

Parameters

```
message the message
```

7.2.2.8 sensordata()

```
bool Motor::sensordata ( )
```

sensordata returns the current PCBMotor optical encoder wheel sendor state

Returns

the current PCBMotor optical encoder wheel sendor state

7.2.2.9 showStatusMessage

showStatusMessage fills the label in the GUI with text

Parameters

age the text to be shown in the label	message
---------------------------------------	---------

7.2.2.10 stop

```
void Motor::stop (
          bool stop ) [slot]
```

stop Tries to stop movenents if possible

Parameters

stop Input: True if movents shall be stopped if possible

7.2.2.11 write

write writes to the serialport

Parameters

data	data to be written to the serial port
	portar in the second portar

7.2.3 Member Data Documentation

7.2.3.1 hometimer

```
QTimer Motor::hometimer [private]
```

Used to iterate through the steps of 'go to the starting position' - but in an inaccurate way.

See also

command_home()

The documentation for this class was generated from the following files:


- motor.h
- · motor.cpp

7.3 UDPlistener Class Reference


The UDPlistener class is used to control dinspect with UDP packages.

```
#include <udplistener.h>
```

Inheritance diagram for UDPlistener:



Collaboration diagram for UDPlistener:



Public Slots

· void bind ()

bind Binds to a new UDP port

Signals

void changeWPangles (QVector< double > angles)

changeWPangles emitted when received message contains command to set waveplate characterization angles ('H')

void changeoffsetusage (bool useoffset)

changeoffsetusage emitted when massage to change usage of offset is received

• void Move (QString controller, double HWPang, double QWPang)

Move emitted when massage to move the waveplates in a cage to certain angles is received.

void MoveHV (QString controller)

MoveHV emitted when massage to move cage to H/V basis is received.

void MovePM (QString controller)

MovePM emitted when massage to move cage to P/M basis is received.

void MoveLR (QString controller)

MoveLR emitted when massage to move cage to R/L basis is received.

Public Member Functions

UDPlistener (QSettings *settings, QObject *parent=0)

UDPlistener listen for commands on a UDP port and execute them.

Private Slots

• void processPendingDatagrams ()

processPendingDatagrams reads data from the UDP socket

• void processCommands (QString msg)

processCommands extracts commands out of received data and executes them

Private Attributes

QSettings * settings

configuration

QUdpSocket socket

the UDP socket

uint port

port the listener listens to

· bool alreadybound

true if the listener is already listening to a port

7.3.1 Detailed Description

The UDPlistener class is used to control dinspect with UDP packages.

Bug There are no known bugs

7.3.2 Constructor & Destructor Documentation

7.3.2.1 UDPlistener()

UDPlistener listen for commands on a UDP port and execute them.

Parameters

settings a pointer to a qsettings instance, used to get the port to bind to and known commands

UDPlistener opens a UDP socket and binds to a port specified in qsettings. Incoming packages are analysed to check if they contain known commands. If they do, these commands are executed.

The commands can be changed in the dinspect settings dialog, but the standard ones are:

- Move(QString, QString)
- Move(QString, double, double)

7.3.3 Member Function Documentation

7.3.3.1 bind

```
void UDPlistener::bind ( ) [slot]
```

bind Binds to a new UDP port

This function checks whether the UDP socket is already bound to a a specific port. If so, it closes this connection and binds to the new port. If not, it binds to the port right away.

7.3.3.2 changeoffsetusage

changeoffsetusage emitted when massage to change usage of offset is received

Parameters

useoffset	true if waveplate offset is to be used
-----------	--

7.3.3.3 changeWPangles

```
void UDPlistener::changeWPangles ( {\tt QVector} < {\tt double} > {\tt angles} \; ) \quad [{\tt signal}]
```

changeWPangles emitted when received message contains command to set waveplate characterization angles ('H')

Parameters

angle	vector containing all angles	. Ordering: HWP0,HWP1,,HWPn,QWP0,QWP1,,QWPn
-------	------------------------------	---

7.3.3.4 Move

Move emitted when massage to move the waveplates in a cage to certain angles is received.

Parameters

controller	either colorcode of cage or 'all'
HWPang	angle of the HWP in degree
QWPang	angle of the QWP in degree

7.3.3.5 MoveHV

MoveHV emitted when massage to move cage to H/V basis is received.

Parameters

controller either colorcode of stage, or 'all'
--

7.3.3.6 MoveLR

```
void UDPlistener::MoveLR (
```

```
QString controller ) [signal]
```

MoveLR emitted when massage to move cage to R/L basis is received.

Parameters

```
controller either colorcode of stage, or 'all'
```

7.3.3.7 MovePM

MovePM emitted when massage to move cage to P/M basis is received.

Parameters

controller	either colorcode of stage, or 'all'
------------	-------------------------------------

7.3.3.8 processCommands

processCommands extracts commands out of received data and executes them

Parameters

msg input: the received message

The documentation for this class was generated from the following files:

- · udplistener.h
- udplistener.cpp

32 Class Documentation


Chapter 8

File Documentation


8.1 debug.h File Reference

contains debug macros

#include "stdio.h"
#include "defines.h"
Include dependency graph for debug.h:



This graph shows which files directly or indirectly include this file:



8.1.1 Detailed Description

contains debug macros


Bug Printing to console does not work on Windows. Workaround: Redirect stderr to stdout and redirect stdout to a file.

This file defines macros to style and simplify output to console.

8.2 defines.h File Reference

Various compile-time definitions.

This graph shows which files directly or indirectly include this file:



Macros

- #define **DEBUGSPECTROMETERCONFIG** FALSE
- #define UNUSED(expr) do { (void)(expr); } while (0)
- #define DEBUG 1
- #define DEBUGERROR 1
- #define DEBUGWARNING 1
- #define DEBUGINFO 1
- #define EPS 0.0000001
- #define PI 3.14159265358979323846
- #define DEGTORAD PI/180
- #define RADTODEG 180/PI

8.2.1 Detailed Description

Various compile-time definitions.

Bug There are no known bugs.

Contains definitions of various kind - mathematical, version constants, debug-variables, \dots

8.2.2 Macro Definition Documentation

8.2.2.1 **DEBUG**

#define DEBUG 1

Enables the execution of various debug-paths used during development.

default: FALSE.

8.2.2.2 DEBUGERROR

#define DEBUGERROR 1

If set to TRUE, enables the execution of the DEBUG_ERROR() maken which is used to write error messages (critical) to stdout.

default: true

8.2.2.3 DEBUGINFO

#define DEBUGINFO 1

If set to TRUE, enables the execution of the DEBUG_INFO() makro which is used to write usefull information to stdout.

default: true

8.2.2.4 DEBUGWARNING

#define DEBUGWARNING 1

If set to TRUE, enables the execution of the DEBUG_WARNING() makro which is used to write warnings about unexpected behaviour to stdout.

default: true

8.2.2.5 DEGTORAD

#define DEGTORAD PI/180

Conversion factor from degree to radians. PI/180

8.2.2.6 EPS

#define EPS 0.0000001

'epsilon' used to check floatingpoint variables in if-conditions.

default: 0.0000001

8.2.2.7 PI

```
#define PI 3.14159265358979323846
```

Pi.

default: 3.14159265358979323846

8.2.2.8 RADTODEG

```
#define RADTODEG 180/PI
```


Conversion factor from radians to degree. 180/PI

8.2.2.9 UNUSED


Use UNUSED(var) in function f(...,type var,...) to silence compiler warnings about unused parameter var.

8.3 helper.h File Reference

```
#include "debug.h"
#include "defines.h"
#include <QMessageBox>
Include dependency graph for helper.h:
```



This graph shows which files directly or indirectly include this file:



Namespaces

helper

contains small functions to display messages

Functions

• void helper::message (QString msg)

message displays a message box

• void helper::error (QString msg)

error displays an error-messagebox and writes a debug_error message to stdout

void helper::warning (QString msg)

warning displays warning-messagebox and writes a debug_warning message to stdout


• void helper::info (QString msg)

info displays an info-messagebox and writes a debug_info message to stdout


8.4 motor.h File Reference

```
#include <QDialog>
#include <QtCore/QtGlobal>
#include <QtSerialPort/QSerialPort>
#include <QObject>
#include <QTimer>
```

#include "defines.h"
Include dependency graph for motor.h:



This graph shows which files directly or indirectly include this file:



Classes


· class Motor

The Motor class operates the PCB-motor.


8.5 udplistener.h File Reference

```
#include <QObject>
#include <QHostAddress>
#include <QSettings>
#include <QUdpSocket>
#include "helper.h"
```

Include dependency graph for udplistener.h:



This graph shows which files directly or indirectly include this file:



Classes

· class UDPlistener

The UDPlistener class is used to control dinspect with UDP packages.

8.6 version.h File Reference

This file contains information about the code version.

Macros

#define VERSION_GIT "v0.1-9-g0c78d25"

The git commit description.

• #define VERSION_GIT_DATE 201812181743

The date of the git commit.

• #define VERSION_BUILD_DATE 201812181755

The builddate.

8.6 version.h File Reference 41

8.6.1 Detailed Description

This file contains information about the code version.

Author

Peter Schiansky

Bug There are no known bugs.

The definitions in this file are used to fill the about-dialog with information about the code: The git commit description, the commit date and the builddate.

Index

bind	PI, 36
UDPlistener, 29	RADTODEG, 37
	UNUSED, 37
cagecontrol, 13	0.10022, 0.
initconnections, 16	EPS
LoadConfig, 17	defines.h, 36
motorGB, 17	error
movemotor, 17	helper, 11
SaveConfig, 17	•
slot_changeWPangles, 18	handleError
slot_changeoffsetusage, 18	Motor, 25
slot moveHV, 18	helper, 11
slot_moveLR, 19	error, 11
	info, 12
slot_movePM, 19	message, 12
slot_movemotors, 19	warning, 12
updatesettings, 19	helper.h, 37
updatesettingsint, 20	hometimer
updatestatus, 20	Motor, 26
useoffset, 20	Wotor, 20
changeWPangles	info
UDPlistener, 30	helper, 12
changeoffsetusage	initconnections
UDPlistener, 29	cagecontrol, 16
command_microstep	isopen
Motor, 23	•
command_moveboth	Motor, 25
Motor, 24	LoadConfig
command_singlestep	cagecontrol, 17
Motor, 24	Cageconiioi, 17
command_step	message
Motor, 24	helper, 12
Motor, 24	Motor, 21
DEBUGERROR	command_microstep, 23
defines.h, 36	command_moveboth, 24
DEBUGINFO	command_singlestep, 24
defines.h, 36	
DEBUGWARNING	command_step, 24 handleError, 25
defines.h, 36	hometimer, 26
DEBUG	isopen, 25
defines.h, 36	motorstatusmessage, 25
DEGTORAD	sensordata, 25
defines.h, 36	showStatusMessage, 26
debug.h, 33	stop, 26
defines.h, 34	write, 26
DEBUGERROR, 36	motor.h, 38
DEBUGINFO, 36	motorGB
DEBUGWARNING, 36	cagecontrol, 17
DEBUG, 36	motorstatusmessage
DEGTORAD, 36	Motor, 25
FPS 36	Move

44 INDEX

UDPlistener, 30	cagecontrol, 20
MoveHV	version.h, 40
UDPlistener, 30 MoveLR	version.n, 40
UDPlistener, 30	warning
MovePM	helper, 12
UDPlistener, 31	write
movemotor	Motor, 26
cagecontrol, 17	
PI	
defines.h, 36	
processCommands	
UDPlistener, 31	
RADTODEG	
defines.h, 37	
SaveConfig	
cagecontrol, 17	
sensordata	
Motor, 25	
showStatusMessage	
Motor, 26	
slot_changeWPangles	
cagecontrol, 18	
slot_changeoffsetusage	
cagecontrol, 18 slot moveHV	
cagecontrol, 18	
slot moveLR	
cagecontrol, 19	
slot movePM	
cagecontrol, 19	
slot_movemotors	
cagecontrol, 19	
stop	
Motor, 26	
UDPlistener, 27	
bind, 29	
changeWPangles, 30	
changeoffsetusage, 29	
Move, 30	
MoveHV, 30	
MoveLR, 30	
MovePM, 31	
processCommands, 31	
UDPlistener, 29	
UNUSED	
defines.h, 37	
udplistener.h, 39	
updatesettings cagecontrol, 19	
updatesettingsint	
cagecontrol, 20	
updatestatus	
cagecontrol, 20	
useoffset	